

Anti-Human/Mouse/Rat FOXG1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11274

3th Edition

Description: This locus encodes a member of the forked-head transcription factor family. The encoded protein, which functions as a repressor, may play a role in brain development. Mutations at this locus have been associated with Rett syndrome.

Antigen: Synthetic peptide of human FOXG1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expression is restricted to the neurons of the developing telencephalon.