

Anti-Human/Mouse AGPAT6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11317

3th Edition

Description: Lysophosphatidic acid acyltransferases (EC 2.3.1.51) catalyze the conversion of lysophosphatidic acid (LPA) to phosphatidic acid (PA). LPA and PA are involved in signal transduction and lipid biosynthesis.

Antigen: Synthetic peptide of human AGPAT6

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitous. Relatively high level of expression in skeletal muscle, heart and testis. Relatively low level of expression in lung.