

**Anti-Human PRSS50 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11320**

---

3th Edition

**Description:** May be involved in proteolysis through its threonine endopeptidase activity. Testis specific. Differentially expressed in some breast cancer tissues.

**Antigen:** Synthetic peptide of human PRSS50

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Testis specific. Differentially expressed in some breast cancer tissues.