

Anti-Human CD66e Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11440

3th Edition

Description:Members of the CEACAM subfamily, including CEACAM5, belong to the CEA gene family. For general information on the CEA gene family, see CEACAM1. Cell surface glycoprotein that plays a role in cell adhesion and in intracellular signaling. Receptor for E.coli Dr adhesins. Found in adenocarcinomas of endodermally derived digestive system epithelium and fetal colon.

Antigen:Synthetic peptide of human CEACAM5

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Found in adenocarcinomas of endodermally derived digestive system epithelium and fetal colon.