

Anti-Human BMP3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11507

3th Edition

Description:BMP3 belongs to the transforming growth factor-beta (TGFB) superfamily. Bone morphogenic protein, also known as osteogenin, induces bone formation.

Antigen:Synthetic peptide of human BMP3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in adult and fetal cartilage.