

Anti-Human/Mouse/Rat RAB3C Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11536

3th Edition

Description: Ras-related protein Rab-3C is a protein that in humans is encoded by the RAB3C gene. RAB proteins, such as RAB3C, are small GTPases implicated in the regulation of vesicular trafficking between membrane compartments. It is involved in protein transport and vesicular traffic probably.

Antigen: Synthetic peptide of human RAB3c

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH 7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1 µg/ml.

Specificity: Expressed in brain, placenta and lung.