

**Anti-Human/Mouse/Rat ANGPT2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11555**

---

3th Edition

**Description:**The protein encoded by this gene is an antagonist of angiotensin 1 (ANGPT1) and endothelial TEK tyrosine kinase (TIE-2, TEK). The encoded protein disrupts the vascular remodeling ability of ANGPT1 and may induce endothelial cell apoptosis. Three transcript variants encoding three different isoforms have been found for this gene.

**Antigen:**Recombinant protein of human ANGPT2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**