

**Anti-Human/Mouse/Rat RAB26 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11571**

---

3th Edition

**Description:** This gene encodes a small GTP-binding protein that belongs to the largest family within the Ras superfamily. RAB26 are important regulators of vesicular fusion and trafficking. The RAB family of small G proteins regulates intercellular vesicle trafficking, including exocytosis, endocytosis, and recycling.

**Antigen:** Synthetic peptide of human RAB26

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Predominantly expressed in brain.