

Anti-Human/Mouse ADIPOQ Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11590

3th Edition

Description: This gene is expressed in adipose tissue exclusively. It encodes a protein with similarity to collagens X and VIII and complement factor C1q. The encoded protein circulates in the plasma and is involved with Human, Mouse metabolic and hormonal processes. Mutations in this gene are associated with adiponectin deficiency. Multiple alternatively spliced variants, encoding the same protein, have been identified.

Antigen: Synthetic peptide of human ADIPOQ

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Synthesized exclusively by adipocytes and secreted into plasma.