

Anti-Human/Mouse STXBP4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11632

3th Edition

Description:The STXBP4 antibody from Proteintech is a rabbit polyclonal antibody to a recombinant protein of human STXBP4. This antibody recognizes human, mouse, rat antigen. The STXBP4 antibody has been validated for the following applications: ELISA, WB, IHC analysis.

Antigen:Recombinant protein of human STXBP4

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: