

Anti-Human/Mouse/Rat ST18 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11666

3th Edition

Description: ST18 (suppression of tumorigenicity 18), also known as ZNF387 ((zinc finger protein 387), is a 1,047 amino acid nuclear protein that is expressed at low levels in heart, liver, kidney, skeletal muscle, pancreas, testis, ovary and prostate. Containing six C2HC-type zinc fingers and an SMC domain, ST18 is a repressor protein that binds to DNA sequences containing a specific bipartite element. ST18 is considered a transcription factor that represses basal transcription activity from target promoters and Inhibits colony formation in cultured breast cancer cells.

Antigen: Recombinant protein of human ST18

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Detected at low levels in heart, liver, kidney, skeletal muscle, pancreas, testis, ovary and prostate. Detected at even lower levels in mammary epithelial cells and breast cancer cells.