

Anti-Human SSBP4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11676

3th Edition

Description:SSBP4 (single stranded DNA binding protein 4) is a 385 amino acid nuclear protein that contains one LisH domain and multiple phosphorylated amino acid residues. The gene encoding SSBP4 maps to human chromosome 19, which consists of over 63 million bases, houses approximately 1,400 genes and is recognized for having the greatest gene density of the human chromosomes. It is the genetic home for a number of immunoglobulin (Ig) superfamily members, including the killer cell and leukocyte Ig-like receptors, a number of ICAMs, the CEACAM and PSG family and Fc receptors (FcRs).

Antigen:Recombinant protein of human SSBP4

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: