

Anti-Human/Mouse SMYD3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11704

3th Edition

Description: This gene encodes a histone methyltransferase which functions in RNA polymerase II complexes by an interaction with a specific RNA helicase. Multiple transcript variants encoding different isoforms have been found for this gene.

Antigen: Recombinant protein of human SMYD3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in skeletal muscles and testis. Overexpressed in a majority of colorectal and hepatocellular carcinomas.