

**Anti-Human/Mouse/Rat SENP8 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11729**

---

3th Edition

**Description:** This gene encodes a cysteine protease that is a member of the sentrin-specific protease family. The encoded protein is involved in processing and deconjugation of the ubiquitin-like protein termed, neural precursor cell expressed developmentally downregulated 8. Alternate splicing results in multiple transcript variants.

**Antigen:** Recombinant protein of human SENP8

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Broadly expressed, with highest levels in kidney and pancreas.