

Anti-Human/Mouse/Rat RBBP9 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11781

3th Edition

Description:The protein encoded by this gene is a retinoblastoma binding protein that may play a role in the regulation of cell proliferation and differentiation. Two alternatively spliced transcript variants of this gene with identical predicted protein products have been reported, one of which is a nonsense-mediated decay candidate.

Antigen:Recombinant protein of human RBBP9

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Widely expressed. Expressed at higher levels in tumor tissues than in normal tissues.