

**Anti-Human/Mouse/Rat PPP2R3C Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11807**

---

3th Edition

**Description:** Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit gamma is an enzyme that in humans is encoded by the PPP2R3C gene. May regulate MCM3AP phosphorylation through phosphatase recruitment. May play a role in the activation-induced cell death of B-cells By similarity.

**Antigen:** Recombinant protein of human PPP2R3C

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitously expressed in brain and other tissues.