

Anti-Human/Mouse PLTP Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11813

3th Edition

Description:The protein encoded by this gene is one of at least two lipid transfer proteins found in human plasma. The encoded protein transfers phospholipids from triglyceride-rich lipoproteins to high density lipoprotein (HDL). In addition to regulating the size of HDL particles, this protein may be involved in cholesterol metabolism. At least two transcript variants encoding different isoforms have been found for this gene.

Antigen:Recombinant protein of human PLTP

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Wide tissue distribution. Placenta > pancreas > lung > kidney > heart > liver > skeletal muscle > brain.