

Anti-Human/Mouse OS9 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11840

3th Edition

Description: This gene encodes a protein that is highly expressed in osteosarcomas. This protein binds to the hypoxia-inducible factor 1 (HIF-1), a key regulator of the hypoxic response and angiogenesis, and promotes the degradation of one of its subunits. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

Antigen: Recombinant protein of human OS9

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitously expressed. Found as well in all tumor cell lines analyzed, amplified in sarcomas. Highly expressed in osteosarcoma OSA-CL and rhabdomyosarcoma RH30 cell lines. Isoform 2 is the major isoform detected in all cell types examined.