

**Anti-Human NUSAP1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11844**

---

3th Edition

**Description:**NUSAP1 is a nucleolar-spindle-associated protein that plays a role in spindle microtubule organization. May associate with chromosomes and promote the organization of mitotic spindle microtubules around them.

**Antigen:**Recombinant protein of human NUSAP1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**