

Anti-Human/Mouse/Rat NDEL1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11847

3th Edition

Description: This gene encodes a coiled-coil protein that plays a role in multiple processes including cytoskeletal organization, cell signaling and neuron migration, outgrowth and maintenance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome X.

Antigen: Recombinant protein of human NDEL1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in brain, heart, kidney, liver, lung, pancreas, placenta and skeletal muscle.