

Anti-Human/Mouse NRIP3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11851

3th Edition

Description:NRIP3 (nuclear receptor interacting protein 3), also known as C11orf14 or NY-SAR-105, is a 241 amino acid protein that is encoded by a gene which Human, Mouse maps to human chromosome 11. With approximately 135 million base pairs and 1,400 genes, chromosome 11 comprises approximately 4% of human genomic DNA and is considered a gene and disease association dense chromosome. The chromosome 11 encoded Atm gene is important for regulation of cell cycle arrest and apoptosis following double strand DNA breaks.

Antigen:Recombinant protein of human NRIP3

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: