

Anti-Human/Mouse/Rat NIT1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11856

3th Edition

Description: This gene encodes a member of the nitrilase protein family with homology to bacterial and plant nitrilases, enzymes that cleave nitriles and organic amides to the corresponding carboxylic acids plus ammonia. Multiple transcript variants encoding different isoforms have been found for this gene.

Antigen: Recombinant protein of human NIT1

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Detected in heart, brain, placenta, liver, skeletal muscle, kidney and pancreas.