

Anti-Human/Mouse CMC4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11891

3th Edition

Description: This gene was identified by involvement in some t(X;14) translocations associated with Human, Mouseature T-cell proliferations. This region has a complex gene structure, with a common promoter and 5' exon spliced to two different sets of 3' exons that encode two different proteins. This gene represents the downstream 8 kDa protein that localizes to mitochondria.

Antigen: Recombinant protein of human CMC4

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in many tissues with a relatively high level in skeletal muscle.