

Anti-Human/Mouse ITK Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11975

3th Edition

Description: This gene encodes an intracellular tyrosine kinase expressed in T-cells. The protein contains both SH2 and SH3 domains which are often found in intracellular kinases. It is thought to play a role in T-cell proliferation and differentiation.

Antigen: Recombinant protein of human ITK

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: T-cell lines and natural killer cell lines.