

Anti-Human/Mouse/Rat FBP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12102

3th Edition

Description: This gene encodes a gluconeogenesis regulatory enzyme which catalyzes the hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate and inorganic phosphate.

Antigen: Recombinant protein of human FBP2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in skeletal muscle (at protein level).