

Anti-Human ARHGEF5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12920

3th Edition

Description:Rho GTPases play a fundamental role in numerous cellular processes initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form a complex with G proteins and stimulate Rho-dependent signals. This protein may be involved in the control of cytoskeletal organization.

Antigen:Recombinant protein of human ARHGEF5

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Mainly expressed in kidney, liver, pancreas, lung and placenta.