

Interferon regulatory factor-3,1-112 aa, Human, Recombinant, E.coli

Cat.NO.: TP02670

3th Edition

Synonyms:IRF3

Description:Members of the Interferon regulatory factor (IRF) family regulate gene expression critical to immune response, hemopoiesis, and proliferation. IRF-3 is a member of the IRF family, and is distinct from other family members. Its transcriptional activity is regulated solely by posttranslational modifications. It plays a crucial role in activation of innate immunity and inflammation in response to viral infection.

Form:Liquid. In Phosphate-buffered saline (pH 7.4)

Molecular Weight:13.0kDa (112aa), confirmed by MALDI-TOF

Sequences:

Purity:> 95% by HPLC

Concentration:1 mg/ml (determined by Bradford assay)

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.