

MRFAP1L1, 1-127aa Human, His tag, E.coli

Cat.NO.: TP03017

3th Edition

Synonyms:MORF4 family-associated protein 1-like, Morf4 family associated protein 1-like 1, PP784

Description:MRFAP1L1 belongs to the MORF4 family-associated protein family. Recombinant human MRFAP1L1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

Form:Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 20% glycerol, 1mM DTT

Molecular Weight:17.2 kDa (150aa) confirmed by MALDI-TOF

Sequences:

MGSSHHHHHHSSGLVPRGSHMGSMRPLDIDEVEAPEEVEVLEPEEDFEQFLLPVINEMREDIASLIREHGRAYLRT
RSKLEWMDNMLIQIKTQVEASEESALNHVQHPSGEADERVSELCEKAEKAKEIAKMAEMLVELVWRIERSESS

Purity:> 95% by HPLC

Concentration:0.5 mg/ml (determined by Bradford assay)

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.