

**NME1, 1-152 aa, Human, Recombinant, E.coli**

**Cat.NO.: TP03165**

---

3th Edition

**Synonyms:**Non-metastatic cells 1, Nucleoside diphosphate kinase A, NDP kinase A, AWD, GAAD, NB, NBS, NDPK-A, NM23, NM23-H1

**Description:**Non-metastatic cells 1 (NME1), also known as NM23-H1, originally identified as a candidate metastasis suppressor gene. NME1 is expressed in different tumor types where their levels have been alternatively associated with reduced or increased metastatic potential. Reductions in NME1 expression have been significantly associated with aggressive behavior in melanoma, breast, colon, and gastric carcinomas.

**Form:**Liquid. In 20mM Tris-HCl buffer (pH7.5) containing 1mM DTT, 10% glycerol

**Molecular Weight:**17.1 kDa (152aa)

**Sequences:**

MANCERTFIAIKPDGVQRGLVGEIIRFEQKGFRLVGLKFMQASEDLLKEHYVDLKDPPFFAGLVKYMHS GPVVAM  
VWEGLNVVKTGRV MLGETNPADSKPGTIRGDFCIQVGRNIIHGSDSVESAEKEIGLWFHPEELVDYTSCAQNWIYE

**Purity:**> 95% by HPLC

**Concentration:**1 mg/ml (determined by Bradford assay)

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.