

**OXLD1, 46-147aa, Human, His tag, E.coli**

**Cat.NO.: TP03273**

---

3th Edition

**Synonyms:**Oxidoreductase-like domain-containing protein 1, C17orf90

**Description:**OXLD1 contains 1 Oxidoreductase-like domain. Recombinant human OXLD1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

**Form:**Liquid. In 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 50% glycerol, 1mM DTT

**Molecular Weight:**13.4 kDa (125aa) confirmed by MALDI-TOF

**Sequences:**

MGSSHHHHHHSSGLVPRGSHMGSPGAQAPDGRRKFGTDHVEVGSQAGADGTRPPKASLPPELQPPTNCCMSG  
CPNCVWVEYADRLLQHFQDGGERALAAL EEHVADENLKAFLRMEIRLHTRCGG

**Purity:**> 95% by HPLC

**Concentration:**0.25 mg/ml (determined by Bradford assay)

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.