

Ppif, 30-206aa, Rat, His tag, E.coli

Cat.NO.: TP03490

3th Edition

Synonyms:Peptidyl-prolyl cis-trans isomerase F, mitochondrial, CyP-D, CypD, PPIase

Description:Ppif also known as Peptidyl-prolyl cis-trans isomerase F, mitochondrial. Ppif accelerates the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and involved in regulation of the mitochondrial permeability transition pore (mPTP). It is proposed that its association with the mPTP is masking a binding site for inhibiting inorganic phosphate (Pi) and promotes the open probability of the mPTP leading to apoptosis or necrosis the requirement of the PPIase activity for this function is debated. Recombinant rat Ppif was expressed in E.coli and purified by using conventional chromatography techniques.

Form:Liquid. In Phosphate buffered saline (pH7.4) containing 10% glycerol, 1mM DTT

Molecular Weight:21.2Da (200aa) confirmed by MALDI-TOF

Sequences:

MGSSHHHHHSSGLVPRGSHMGSCSDGGARGANSSSQNPLVYLDVGADGQPLGRVVLELKADVVPKTAENFRA
LCTGEKGFYKYSTFHRVIPAFMCQAGDFTNHNGTGGKSIYGSRFPDENFTLKHVGPVLSMANAGPNTNGSQFF
ICTIKTDWLDGKHVVFGHVKEGMDVVKKIESFGSKSGKTSKKIVITDCGQLS

Purity:> 95% by HPLC

Concentration:1mg/ml (determined by Absorbance at 280nm)

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.