

Protein Marker-1

Cat.NO.: TP03541

3th Edition

Synonyms:n/a

Description:Mid-Range Protein Marker-1 consists of Seven precisely sized recombinant proteins (M.W. 70, 57, 40, 28, 18, 13.5, 8.5KDa). The concentration of each protein is 0.1ug/ul except for 8.5KDa (0.2ug/ul). Usage Note: Use 5ul per lane and you need not heat before use.

Form:12mM Tris-HCl pH6.8, 5% glycerol, 0.4% SDS 2.88 mM 2-mercaptoethanol, 0.02% bromophenol blue

Molecular Weight:34

Sequences:

Purity:> 95% by HPLC

Concentration:n/a

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.