

Protein Marker-2, Broad range

Cat.NO.: TP03542

3th Edition

Synonyms:

Description: ATGen Broad-Range Protein Marker-2 consists of eight precisely sized recombinant proteins (M.W. 150, 100, 70, 50, 35, 25, 20, 15 KDa). The concentration of each protein is about 0.1 µg/µl. For research use only, not for use in diagnostic procedures.

Form: Please see datasheet.

Molecular Weight:

Sequences:

Purity: > 95% by HPLC

Concentration: Please see datasheet.

Endotoxin Level: <1.0 EU per 1 µg of protein (determined by LAL method)

Storage: Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.