

UBE2J2, 1-226aa, Human, His tag, E.coli

Cat.NO.: TP04388

3th Edition

Synonyms:ubiquitin-conjugating enzyme E2 J2, NCUBE-2, NCUBE2, PRO2121

Description:UBE2J2, also known as ubiquitin-conjugating enzyme E2 J2, belongs to the ubiquitin-conjugating enzyme family. The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is located in the membrane of the endoplasmic reticulum. Recombinant human UBE2J2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.

Form:Liquid. In 20 mM Tris-HCl buffer, pH8.0, 30% glycerol, 1mM DTT, 200mM NaCl

Molecular Weight:28.0kDa (250aa) confirmed by MALDI-TOF

Sequences:

MGSSHHHHHSSGLVPRGSHMGSHMSSTSSKRAPTTATQRLKQDYLRICKDPVPYICAEPLPSNILEWHYVVRGP
EMTPYEGGYH GKLI FPREFPFKPPSIYMITPNGRFKCNTRLCLSITDFHPDTWNPAWSVSTILTGLLSFMVEKGPTL
GSIETSDFTKRQLAVQSLAFNLKDKVFCELFPEVVEEIKQKQKAQDELSSRPQTLPLPDVVPDGETHLVQNGIQLLN
GHAPGAVPNLAGLQQANRHH

Purity:> 95% by HPLC

Concentration:1mg/ml (determined by Bradford assay)

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.