

**Recombinant Human Neuritin Protein**

**Cat.NO.: TP05089**

3th Edition

**Synonyms:**Neuritin;NRN1;NRN

**Description:**Neuritin/NRN1 is a member of the neuritin family and can be expressed in postmitotic-differentiating neurons of the developmental nervous system and neuronal structures associated with plasticity in the adult. Neuritin/NRN1 promotes neurite outgrowth, arborization and neuritogenesis. The protein contains a consensus cleavage signal found in glycosylphosphatidylinositol (GPI)-anchored proteins. Overexpression of the encoded protein may be associated with astrocytoma progression.

**Form:**PBS

**Molecular Weight:**12.1 kDa

**Sequences:**Ala28-Gly116

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.