

Recombinant Human ADAM-like Protein Decysin-1/ADAMDEC1 Protein(C-6His)

Cat.NO.: TP05093

3th Edition

Synonyms:ADAM DEC1;a disintegrin and metalloproteinase domain-like protein decysin-1;ADAM-like protein decysin-1

Description:ADAM DEC1 protein is expressed highly in the small intestine and appendix, moderately in lymph node, mucosal lining of the colon, thymus, spleen and very weakly in the bone marrow. ADAM DEC1 is induced during DC maturation and up-regulated in response to T-cell signals. It may play an important role in the control of the immune response and during pregnancy.

Form:PBS

Molecular Weight:30.3 kDa

Sequences:Lys206-Glu470

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.