

Recombinant Human Calreticulin-3/CALR3/CRT2 Protein

Cat.NO.: TP05185

3th Edition

Synonyms:Calreticulin-3; calreticulin-2; calsperin; CALR3; CRT2

Description:Calreticulin-3 belongs to the calreticulin family, members of which are calcium binding chaperones localized mainly in the endoplasmic reticulum. It can be divided into a N-terminal globular domain, a proline-rich P-domain forming an elongated arm-like structure and a C-terminal acidic domain. During spermatogenesis process, Calreticulin-3 may act as a lectin-independent chaperone for specific client proteins such as ADAM3. Defects in CALR3 are the cause of familial hypertrophic cardiomyopathy type 19 (CMH19), it is a hereditary heart disorder characterized by ventricular hypertrophy, which is usually asymmetric and often involves the interventricular septum. The symptoms include dyspnea, syncope, collapse, palpitations, and chest pain.

Form:PBS

Molecular Weight:42.9 kDa

Sequences:Thr 20-Leu384

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.