

Recombinant Human Nucleolar Protein 3/NOL3 Protein(N-GST)

Cat.NO.: TP05278

3th Edition

Synonyms:Nucleolar Protein 3; Apoptosis Repressor With CARD; Muscle-Enriched Cytoplasmic Protein; Myp; Nucleolar Protein of 30 kDa; Nop30; NOL3; ARC; NOP

Description:Nucleolar Protein 3 is encoded by NOL3 gene; multiple transcript variants encoding different isoforms have been found for this gene. So far, Nucleolar protein 3 has show to have two Isoforms. Isoform 1 may be involved in RNA splicing.Isoform 2 may inhibit apoptosis.It has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53.

Form:PBS

Molecular Weight:48.9 kDa

Sequences:Met 1-Ser208

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.