

**Recombinant Human Host Cell Factor 2/HCFC2 Protein(N-6His)**

**Cat.NO.: TP05280**

3th Edition

**Synonyms:**Host cell factor 2 ;HCF-2;C2 factor;Hcfc2

**Description:**Host cell factor 2(HCFC2) is a cytoplasmic protein. It contains 2 fibronectin type-III domains.HCFC2 binds KMT2A/MLL1, as component of the MLL1/MLL complex.Hcfc2 negative regulation of transcription from RNA polymerase II promoter.

**Form:**PBS

**Molecular Weight:**55.2 kDa

**Sequences:**Met 1-Glu497

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.