



Recombinant Human 4-1BB Ligand/4-1BBL/TNFSF9/CD137L Protein(C-6His)

Cat.NO.: TP05336

3th Edition

Synonyms:Tumor necrosis factor ligand superfamily member 9; 4-1BB ligand; 4-1BBL; TNFSF9

Description:Tumor necrosis factor ligand superfamily member 9(4-1BBL) is single-pass type II membrane protein which is a member of the the tumor necrosis factor family. 4-1BBL is a 254 amino acids cytokine that is expressed in brain, placenta, lung, skeletal muscle and kidney. TNFSF9 has been shown to reactivate anergic T lymphocytes in addition to promoting T lymphocyte proliferation. This cytokine may have a role in activation-induced cell death (AICD) and cognate interactions between T-cells and B-cells/macrophages. It has also been shown to be required for the optimal CD8 responses in CD8 T cells, and is thought to be involved in T cell-tumor cell interaction.

Form:PBS

Molecular Weight:41.8 kDa

Sequences:Arg71-Glu254

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.