



**Recombinant Human Protein Disulfide-Isomerase A4/PDIA4 Protein(C-6His)**

**Cat.NO.: TP05434**

3th Edition

**Synonyms:**Protein Disulfide-Isomerase A4; Endoplasmic Reticulum Resident Protein 70; ER Protein 70; ERp70; Endoplasmic Reticulum Resident Protein 72; ER Protein 72; ERp-72; ERp72; PDIA4; ERP70; ERP72

**Description:**Protein Disulfide-Isomerase A4 (PDIA4) is an endoplasmic reticulum luminal protein that belongs to the protein disulfide isomerase family. Human PDIA4 is synthesized as a 625 amino acid precursor that contains a 20 amino acid signal sequence, and a 625 amino acid mature chain, including three thioredoxin domains. PDIA4 catalyzes the rearrangement of -S-S- bonds in proteins and is thought to be a deoxycytidine kinase. In addition, PDIA4 serves as a proteases protein disulfide isomerase, phospholipase or an arrangement of these.

**Form:**PBS

**Molecular Weight:**71.7 kDa

**Sequences:**Val21-Leu645

**Purity:**> 95% by HPLC

**Concentration:**

**Endotoxin Level:**<1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:**Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.