
Recombinant Human LYVE-1/HAR/XLKD1 Protein(C-6His)**Cat.NO.: TP05505**

3th Edition

Synonyms:Lymphatic Vessel Endothelial Hyaluronic Acid Receptor 1; LYVE-1; Cell Surface Retention Sequence-Binding Protein 1; CRSBP-1; Extracellular Link Domain-Containing Protein 1; Hyaluronic Acid Receptor; LYVE1; CRSBP1; HAR; XLKD1

Description:Lymphatic Vessel Endothelial Hyaluronic Acid Receptor 1 is a single-pass type I membrane protein. LYVE-1 is a CD44 homolog found primarily on lymphatic endothelial cells. LYVE-1 is mainly expressed in endothelial cells lining lymphatic vessels. While LYVE-1 functions as a Ligand-specific transporter trafficking between intracellular organelles (TGN) and the plasma membrane. LYVE-1 plays a role in autocrine regulation of cell growth mediated by growth regulators containing cell surface retention sequence binding (CRS). It may act as an hyaluronan (HA) transporter, either mediating its uptake for catabolism within lymphatic endothelial cells themselves, or its transport into the lumen of afferent lymphatic vessels for subsequent re-uptake and degradation in lymph nodes.

Form:PBS

Molecular Weight:24.6 kDa

Sequences:Leu20-Thr 238

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.