

Recombinant Human Follistatin-Related Protein 1/FSTL1 Protein(C-6His, E. coli)

Cat.NO.: TP05516

3th Edition

Synonyms:Follistatin-Related Protein 1; Follistatin-Like Protein 1; FSTL1; FRP

Description:Follistatin-Related Protein 1 (FSTL1) is a secreted protein that contains two EF-hand domains, one follistatin-like domain, one Kazal-like domain, and one VWFC domain. Its functional significance in physiological and pathological processes is not completely understood. However, FSTL1 is thought to modulate the action of some growth factors on cell proliferation and differentiation. FSTL1 maybe an autoantigen associated with rheumatoid arthritis.

Form:PBS

Molecular Weight:33.8 kDa

Sequences:Glu21-Ile308

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.