

Recombinant Human TXNRD1 / TRXR1 Protein (aa 161-647, His tag)

Cat.NO.: TP05589

3th Edition

Synonyms:GRIM-12;TR;TR1;TRXR1;TXNR

Description:Thioredoxin reductase 1 (TXNRD1) which is a selenocysteine-containing protein is overexpressed in many malignancies. TXNRD1 plays a key role in regulating cell growth and transformation, and protects cells against oxidative damage. We investigated the association between TXNRD1 polymorphisms and ATDH susceptibility. Moreover, TXNRD1 is an essential selenium-containing enzyme involved in detoxification of reactive oxygen species (ROS) and redox signaling. And genetic variations in TXNRD1 favor the development of Drug-induced liver injury (DILI), which is the most common adverse drug reaction.

Form:PBS

Molecular Weight:55 kDa

Sequences:Tyr 161-Cys 647

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.