

Recombinant Mouse IgG3-Fc Protein

Cat.NO.: TP05606

3th Edition

Synonyms: IgG3

Description:

Form: PBS

Molecular Weight: 26 kDa

Sequences: Glu 97-Lys 329

Purity: > 95% by HPLC

Concentration:

Endotoxin Level: <1.0 EU per 1 ug of protein (determined by LAL method)

Storage: Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.