

**Recombinant Mouse IgG2b-Fc Protein**

**Cat.NO.: TP05607**

3th Edition

**Synonyms:** IgG2b

**Description:**

**Form:** PBS

**Molecular Weight:** 26.3 kDa

**Sequences:** Glu 97-Lys 335

**Purity:** > 95% by HPLC

**Concentration:**

**Endotoxin Level:** <1.0 EU per 1 ug of protein (determined by LAL method)

**Storage:** Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.