

Recombinant Human HDGFRP3/HRP-3/HDGF2 Protein(C-6His)

Cat.NO.: TP05685

3th Edition

Synonyms:Hepatoma-Derived Growth Factor-Related Protein 3; HRP-3; Hepatoma-Derived Growth Factor 2; HDGF-2; HDGFRP3; HDGF2

Description:Hepatoma-Derived Growth Factor-Related Protein 3 (HDGFRP3) belongs to the HDGF family. HDGFRP3 can be found in testis, heart, spinal cord and brain. HDGFRP3 localizes to the nucleus and contains one PWWP domain. HDGFRP3 enhances DNA synthesis and may have a role in cell proliferation.

Form:PBS

Molecular Weight:23.6 kDa

Sequences:Met 1-Thr 203

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.