

Recombinant Human Neurotrimin/NTRI Protein(C-6His)

Cat.NO.: TP05695

3th Edition

Synonyms:Neurotrimin;NTM; IgLON family member 2; IGLON2; HNT;NTRI;

Description:Neurotrimin is a glycoprotein. It belongs to the immunoglobulin superfamily and IgLON family. Neurotrimin localizes to the cell membrane and contains three Ig-like C2-type domains. Neurotrimin acts as a glycosylphosphatidylinositol (GPI)-anchored neural cell adhesion molecule. Neurotrimin may promote neurite outgrowth and adhesion via homophilic and heterophilic interactions.

Form:PBS

Molecular Weight:32.3 kDa

Sequences:Gly34-Leu316

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.