

Recombinant Mouse Interleukin-21/IL-21 Protein

Cat.NO.: TP05843

3th Edition

Synonyms:Interleukin-21;IL21

Description:Interleukin-21 also known as IL-21 is a protein that in mouse is encoded by the IL21 gene, belongs to the IL-15/IL-21 family. Interleukin-21 is a cytokine that has potent regulatory effects on cells of the immune system, including natural killer (NK) cells and cytotoxic T cells that can destroy virally infected or cancerous cells. This cytokine induces cell division/proliferation in its target cells.

Form:PBS

Molecular Weight:14.4 kDa

Sequences:Pro25-Ser146

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.