

Recombinant Human Polypeptide GalNac Transferase 3/GALNT3 Protein(C-6His)

Cat.NO.: TP05907

3th Edition

Synonyms:Polypeptide N-acetylgalactosaminyltransferase 3; Polypeptide GalNAc transferase 3; GalNAc-T3; pp-GaNTase 3; Protein-UDP acetylgalactosaminyltransferase 3; UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 3; HFTC; HHS

Description:Polypeptide N-acetylgalactosaminyltransferase 3(GALNT3) belongs to the glycosyltransferase 2 family and galNAc-T subfamily. It expressed in organs that contain secretory epithelial glands and it highly expressed in pancreas, skin, kidney and testis. There are two conserved domains in the glycosyltransferase region: the N-terminal domain (domain A, also called GT1 motif), which is probably involved in manganese coordination and substrate binding and the C-terminal domain (domain B, also called Gal/GalNAc-T motif), which is probably involved in catalytic reaction and UDP-Gal binding .This protein plays a major role in regulating phosphate levels within the body (phosphate homeostasis). Among its many functions, phosphate plays a critical role in the formation and growth of bones in childhood and helps maintain bone strength in adults.

Form:PBS

Molecular Weight:69.1 kDa

Sequences:Gln38-Asp633

Purity:> 95% by HPLC

Concentration:

Endotoxin Level:<1.0 EU per 1 ug of protein (determined by LAL method)

Storage:Can be stored at +4°C short term (1-2 weeks). For long term storage, aliquot and store at -20°C or -70°C. Avoid repeated freezing and thawing cycles.